

# CI-602

## ■ Narrow Gauge Root Imager

---


A portable minirhizotron system  
compatible with 2-inch diameter tubes

Capture high-resolution, digital images of living roots in soil over multiple growing seasons with the CI-602 Narrow Gauge Root Imager. Designed for use with 2-inch tubes, users now have the ability to work within new or existing tube installations.





For more information or to receive pricing, please contact [sales@cid-inc.com](mailto:sales@cid-inc.com).



## Product Features

- ▶ Compatible with transparent root tubes which have an inner diameter of two inches (including Bartz tubes).

- ▶ Comes with six durable transparent root tubes:

Available in lengths of 105cm or 182cm

Made of extruded acrylic

Feature custom-machined watertight bottom plug and insulated top cap to prevent moisture build-up in the tube

- ▶ Durable and field-ready
- ▶ Lightweight and easy to transport to any location
- ▶ Suitable for sites with root tubes across a range of treatments or conditions
- ▶ Controlled and powered by included tablet pre-loaded with RootSnap Analysis Software. Use RootSnap to quickly calculate parameters from images including total and individual root:

length

area

volume

diameter

branching angle