

www.felixinstruments.com sales@felixinstruments.com +1 (360) 833-8835



F-750 Produce Quality Meter

The F-750 Produce Quality Meter non-destructively and rapidly analyzes internal traits of fresh produce through near infrared analysis (NIR). The F-750 sends particles of light into the fruit and then measures the light interactance with molecules to assess the quantity of key quality traits. The model-building software included with the device allows for customizable calibration for specific varieties and regions.

Measure produce quality quickly and non-destructively!

Product Features

- Measures traits including:
 - Total Soluble Solids (TSS or brix) Dry Matter Titratable Acidity Internal Color
 - External Color
- Light-weight, portable, and field-ready
- Fast and consistent—scans take 4-6 seconds per sample
- Non-destructive measurements
- Built-in GPS for easy crop mapping

- Transflective display for outdoor viewing
- Rechargeable/replaceable batteries that last all day
- Wi-Fi capabilities
- Ability to create unlimited customizable models for specific varieties or locations
- The F-750 supports English, and Portuguese on the instrument display.
- Use the optional Liquids & Powders Kit to make measurements on commodities in the liquid or powder state

F-750 Specifications

Spectrometer	Carl Zeiss MMS-1 Spectrometer
Range	310-1100nm
Spectral sample size	3nm
Spectral resolution	8-13nm
Light source	Xenon Tungsten Lamp
Lens	Glass, coated to enhance NIR
Shutter	Gold-plated reference standard
Display	Sunlight visible transflective LCD screen
PC interface	USB and Wi-Fi SD card
Data recorded with each measurement	Raw Data, Reflectance, Absorbance, First Derivative Absorbance, Second Derivative Absorbance
Power source	Removable 3100 milliamp hour lithium-ion battery
Battery life	1600+ Measurements
Data storage	Removable 16 GB Wi-Fi SD card
Body	Heavy-duty anodized aluminum body
Weight	1.05 kg

Make a Custom Model

Your F-750 comes with pre-programmed demo models for a range of commodities. Try one of these models to familiarize yourself with the F-750, or create a custom model for precise measurements.



PLIED FOOD SCIEN

Copyright 🕲

2016 Felix Instruments



Applications

 Optimization of crop management and harvest-timing for fresh produce growers:

Mangoes, avocados, apples, pears, citrus, tomatoes, stone fruit, blackberries, cherries, persimmons, grapes and many more!

- Postharvest quality management in cold storage and ripening rooms
- Import quality assessment
- Retail outlet inspection

Visit our website to learn more about how the F-750 is being used around the world.

www.felixinstruments.com sales@felixinstruments.com

Phone: +1 (360) 833-8835 Toll Free: 1-800-767-0119 Fax: +1 (360) 833-1914