

Accessories



LC-1
■ Square Leaf Chamber

For open-system measurements
of trees, shrubs and herbs
with small, broad leaves.
25 mm x 25 mm



LC-2
■ Wide Rectangular
Leaf Chamber

For open-system measurements
of trees, shrubs and herbs
with small, regular leaves.
55 mm x 20 mm



LC-3
■ Narrow Rectangular
Leaf Chamber

For open-system measurements
of grasses and grass-like leaves.
65 mm x 10 mm




LC-4
■ Small Cylindrical
Leaf Chamber

For open-system measurements
of succulents and small-
needled conifers.
25 mm x 90 mm




LC-5
■ Large Cylindrical
Leaf Chamber

For open-system measurements
of large-needled conifers.
50 mm x 70 mm



LC-7
■ 1/4 Liter Leaf Chamber


For closed-system measurements
of long, narrow leaves.
104 mm x 33 mm x 73 mm



LC-8

■ **1/2 Liter Leaf Chamber**


For closed-system measurements of medium-sized leaves.
89 mm x 66 mm x 86 mm



LC-9

■ **1 Liter Leaf Chamber**

For closed-system measurements of large leaves.
112 mm x 91 mm x 99 mm



LC-10

■ **4 Liter Leaf Chamber**

For closed-system measurements of very large leaves.
180 mm x 130 mm x 170 mm



LC-11

■ **Cactus Leaf Chamber**

For measuring the leaves of Cacti with the CI-340 Handheld Photosynthesis System.
Window size: 14.6cm²



CI-301SR

■ **Soil Respiration Chamber**

For measurements of CO₂ flux from soils.



CI-301CC

■ **Canopy Chamber Attachment**

Interface between chamber and instrument which allows user to design and build their own custom chamber.



F-9002L

■ **Sample Chamber**

For measuring small fruits. Includes a solid bottom plate for closed air-tight measurement.



F-9008L

■ **Sample Chamber**

Air-tight chamber for measuring large fruits.



F-900CW

■ **Attachment Wand**

Interface for the F-900 and any sample chamber. Also included in F-900RK Research Kit.